7.1 Determination of moisture

Apparatus

Aluminium moisture cup with lid, hot air oven, metal tongs, desiccator, analytical balance.

Procedure

- Dry aluminium cup in oven at 100°C for 15 to 30 minutes.
- Cool in a dessicator, weigh and record its weight.
- Weight accurately about 5 g of previously ground sample in an aluminium cup and record its weight.
- Place the cup with sample in hot air oven (Fig. 7.1) maintained at 105 ± 2°C and dry for at least 2 h.
- Cool in a desiccator and weigh.
- Repeat the process of heating, cooling and weighing until the difference between two successive weighings is less than one mg.



Fig. 7.1 Hot air oven

Calculation

 $100 (M_1 - M_2)$

Moisture, per cent by mass = ------M₁ M

Where,

 M_1 =mass in g of the cup with the material before drying M_2 =mass in g of the cup with the material after drying and

M = mass in g of the empty cup

Reference: IS 2052-2012, Indian Standards Institution, Specification for compounded feeds for cattle, Fifth revision.